

1 IHE Integration Statement

IHE Integration Statement		Date	2007-01-19
Vendor	Product Name	Version	
INFINITT	Infinitt PACS	3.0.3.2	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager	None	
	Image Display	None	
	Evidence Creator	None	
Patient Information Reconciliation(PIR)	Image Manager	None	
	Report Manager	None	
Consistent Presentation of Images (CPI)	Image Manager	None	
	Image Display	None	
	Evidence Creator	None	
	Print Composer	None	
Presentation of Grouped Procedures (PGP)	Image Manager	None	
Access to Radiology Information(ARI)	Image Manager	None	
	Image Display	None	
	Report Reader	None	
	Report Repository	None	
Key Image Note(KIN)	Image Manager	None	
	Evidence Creator	None	
	Image Display	None	
Simple Image and Numeric Report(SINR)	Report Creator	None	
	Report Manager	None	
	Report Reader	None	
	Report Repository	None	
Basic Security(SEC)	Image Manager	None	
	Image Display	None	
	Print Composer	None	
	Report Creator	None	

	Report Reader	None
	Report Manager	None
	Report Repository	None
	Audit Record Repository	None
	Secure Node	None
Post-Processing Workflow(PWF)	Image Display	None
	Evidence Creator	None
Reporting Workflow(RWF)	Report Manager	None
	Report Creator	None
	Report Reader	None
Evidence Document(ED)	Image Manager	None
	Evidence Creator	None
	Image Display	None
Retrieve Information for Display(RID)	Display	None
	Information Source	None
Patient Identifier Cross-referencing for MPI(PIX)	Patient Identity Consumer	None
Patient Synchronized Application(PSA)	Patient Context Participant	None
Consistent Time(CT)	Time Client	None
Audit Trail and Node Authentication (ATNA)	Audit Record Repository	None
	Secure Node	None
Cardiac Catheterization Workflow (CATH)	Image Manager	None
	Image Display	None
Echocardiography Workflow (ECHO)	Image Manager	None
	Image Display	None
Cross-enterprise Document Sharing (XDS)	Document Consumer	None
	Document Source	None
Portable Data for Imaging (PDI)	Image Display	None
	Print Composer	None
	Report Reader	None
	Display	None
	Portable Media Creator	None
	Portable Media Importer	None

Cross-enterprise Document Sharing – Image (XDS-I)	Document Consumer	None
	Imaging Document Consumer	None
	Imaging Document Source	None
<u>Internet address for vendor's IHE information:</u>		
Not Available		
Links to Standards Conformance Statements for the Implementation		
DICOM	Not Available	
Links to general information on IHE		
In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.iira-net.or.jp/ihe-j